



DRUG REFERENCE TABLE

DRUG NAME	Test	COMMON NAMES	SPECIFIC FORM OF SUBSTANCE	TIME DETECTABLE IN URINE
Amphetamines	5-Panel	Uppers, Speed, Meth, Dexies, Crank, Benny, Jelly Beans, Crystal, White Cross, Black Beauties	Adderall, Diet pills, Methamphetamine, Ritalin	1-2 Days
Cocaine	5-Panel	Coke, Crank, Snow, Blow, Nose Candy, Freebase, C, Flake, Tot, Rock	Cocaine	1-3 Days
Marijuana	5-Panel	Dope, Pot, Reefer, Joint, Grass, Mary Jane	Marijuana	Occasional Use: 1-3 days Chronic Use: Up to 30 days
Opiates	5-Panel	Heroin (Smack, Horse, Junk, China White), Morphine (M)	Morphine, Poppy Seeds, Tylenol #3	1-3 Days
Phencyclidine	5-Panel	Angel Dust, Devil, Stick, PCP, Dummy Dust, Wack	Angel Dust	Occasional Use: 1-5 days Chronic Use: Up to 30 days
Barbiturates	7-Panel	Downers, Dolls, Red, Tuinal, Rainbows, Yellows, Blues, Goofballs, M&M's, Red Devils	Butalbital, Phenobarbital	Pentobarbital: 1-3 days Phenobarbital: 1-3 weeks
Benzodiazepines	7-Panel	Downers, Tranqs	Librium, Valium, Xanax	1-14 days
Methadone	10-Panel	Done, Dolophine, Methadose Dolls	Heroin substitute, Drug treatment for withdrawal symptoms	1-3 days
Methaqualone	10-Panel	Ludes, Sopor, Mequin, Lennon, Quaalue	Quaaludes	1-7 days
Propoxyphene	10-Panel	No Known Names	Darvocet, Darvon	1-3 days

Creatinine – The Creatinine level measures the kidneys ability to filter urine. Urine creatinine is used to evaluate sample collections for adulteration. As a result, greater than or equal to .20 GM/L is normal. A result less than .20 GM/L may indicate possible adulteration.

Specific Gravity – Is the weight of a substance compared with the weight of an equal amount of some other substance taken as a standard. For liquids, the usual standard is water. The Specific Gravity for water is 1.000 and urine is .003 – 1.030.

